

Product specification

Cuttec S · Manual sliding microtome

Item No.: 10040001



Specifications

Operating temperature range	+10 to +35 °C
Operating humidity	max. rel. 80 % non-condensing
Storage temperature range	+5 to +55 °C
Storage humidity	max. rel. 80 % non-condensing
Section thickness range	0 - 60 µm
Section thickness settings	from 0 - 15 µm in 1 µm steps from 15 - 60 µm in 5 µm steps
Automatic feed	0 - 60 µm
Specimen feed	30 mm
Maximum specimen size	50 x 70 x 20 mm
Specimen orientation	± 8 °
Length of cutting way	180 mm
Dimensions (W x D x H)	150 x 450 x 200 mm
Weight (w/o accessories)	approx. 18 kg

Consumables

	Item No.
Disposable blades (plasma coated) Low profile plasma LPS for soft specimens, 50 pcs. / box	28407005
Disposable blades (plasma coated) Low profile plasma LPH for hard specimens, 50 pcs. / box	28407004
Disposable blades (plasma coated) High profile plasma HPS for soft specimens, 50 pcs. / box	28407007
Disposable blades (PTFE coated) Low profile suitable for all types of specimens, 50 pcs. / box	28407000
Disposable blades (PTFE coated) High profile suitable for all types of specimens, 50 pcs. / box	28408000

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
M: mail@slee.de



www.slee.de



[YouTube](#)



[LinkedIn](#)



[Instagram](#)