

Product specification

Paratec Premium, Paratec · Tissue embedding paraffin

various



Description

Paratec Premium
Premium paraffin (pure paraffin with added polymer),
melting range 54 - 57 °C

Quantity

Item No.

1 kg 30010004

Paratec Premium
Premium paraffin (pure paraffin with added polymer),
melting range 54 - 57 °C

2.5 kg 30010006

Paratec Premium
Premium paraffin (pure paraffin with added polymer),
melting range 54 - 57 °C

10 kg 30010008

Paratec
Paraffin (without polymer additive),
melting range 54 - 56 °C

20 kg 30010002

Comments / additional information:

In case of further questions, please feel free to contact us at every time.

SLEE medical GmbH · Am Neuberg 14
55268 Nieder-Olm · Germany
www.slee.de

T: +49 (0) 6136 76997-0
E: mail@slee.de



www.slee.de



YouTube



LinkedIn



Instagram